

# LIGAMENT DEFICIENCIES

BRACE	LIGAMENT DEFICIENCIES SUPPORTED					BRACE OPTIONS/CHARACTERISTICS			
	ACL	MCL	PCL	LCL	COMBINED INSTABILITY (ACL & PCL)	HINGE	TYPE	STRAPPING OPTIONS	CALF CUFF OPTIONS
DEFIANCE	○	○	○	○	○	Polycentric FourcePoint Drop-Lock OA	Rigid, customised	ACL, PCL or CI	Standard/Short
FULLFORCE	○					FourcePoint	Rigid	ACL	Standard/Short
RENEGADE	○	○		○	○	Polycentric	Rigid	ACL or CI	Standard/Short
ARMOR	○	○	○	○	○	Polycentric	Rigid	ACL, PCL or CI	Short
ARMOR FOURCEPOINT	○	○		○		FourcePoint	Rigid	ACL	Short
4-TITUDE			○			Polycentric	Rigid	PCL	Standard/Short
PLAYMAKER II	○	○		○		Polycentric	Soft, wraparound or sleeve	N/A	N/A
PLAYMAKER XPERT		○		○		Polycentric	Soft, wraparound	N/A	N/A
DELUXE HINGED KNEE		○		○		Polycentric	Soft, wraparound	N/A	N/A
REDDIE BRACE		○		○		Polycentric	Soft, wraparound	N/A	N/A

## MEASURING GUIDE



WATCH THE MEASURING VIDEO FOR ASSISTANCE

Measure at three points on the leg and use the chart to determine the size required.

Thigh 15cm above mid-patella (kneecap) · Knee, mid-patella (kneecap) · Calf 15cm below mid-patella (kneecap)

If one of the three measurements is a size up, select that size

For example, Thigh Size 2, Knee Size 3, Calf Size 2. Order Size 3

Note that for OA braces, the off-loading hinge is located on the opposite side of the affected compartment.

SIZE	THIGH CIRC. 15CM ABOVE MID-PATELLA	KNEE CENTRE CIRC. MID-PATELLA	CALF CIRC. 15CM BELOW MID-PATELLA
1 (XS)	33 – 39 cm	30.5 – 33cm	25.5 – 30.5cm
2 (S)	39 – 47 cm	33 – 35.5cm	30.5 – 35.5cm
3 (M)	47 – 53 cm	35.5cm – 38cm	35.5cm – 40.5cm
4 (L)	53 – 60 cm	38 – 43cm	40.5 – 47cm
5 (XL)	60 – 67 cm	43 – 48.25cm	57 – 50.75cm
6 (XXL)	67 – 75 cm	48.25cm – 53.25cm	50.75cm – 56cm
7 (XXXL)	75 – 81 cm	53.25 – 58.5cm	56 – 61cm